

Laboratory Accessories

Boston Round Amber Bottles With Black Phenolic Cap & Inert Liner

<u>Part Number</u>	<u>Bottle Size</u>	<u>Per Case</u>
DCG601A015242	1 oz.-mL	24
DCG602A015242	2 oz.-mL	24
DCG604A015242	4 oz.-mL	24
DCG608A015242	8 oz.-mL	12
DCG616A015242	16 oz.-ml	12
DCGLPLAB032100001	32 oz.-mL	12

Gas Chromatography Vials, Caps & Septa

For GC autosampler or manual injection of samples with Inert septa that will not react with sample

<u>Part Number</u>	<u>Vial Size</u>	<u>Per Case</u>
DCG091200	2 mL-Clear	100
DCG092200	2 ML-Amber	100
DCG90030B	<u>CAPS & SEPTA</u>	100

Gas Chromatography Syringes & Replacement Needles

For manual injection

<u>Part Number</u>	<u>Syringe Size</u>	<u>Per Pack</u>
DCG002030	10µ fixed needle	1
DCG003000	25µ fixed needle	1
DCG003200	25µ gastight	1
DCG002050	10µ removable needle	1
DCG003050	25µ removable needle	1
DCG037110	10µ replacement needle	5 pk
DCG038110	25µ replacement needle	5 pk

DCG offers laboratory accessories for your laboratory needs. These accessories are used in the day to day operation of the laboratory and are necessary for successful sample preparation and analysis.

Nitrile Gloves

For safety precautions during sample prep or handling chemicals Latex free

<u>Part Number</u>	<u>Size</u>	<u>Per Case</u>
DCG19-130-1597A	Small	1000
DCG190-130-1597B	Medium	1000
DCG190-130-1597C	Large	1000
DCG190-130-1597D	Extra Large	100

Disposable Transfer Pipettes

For transferring samples or reagents into bottles or vials made of inert low density polyethylene

<u>Part Number</u>	<u>Size</u>	<u>Per Case</u>
DCG263	Small	100
DCG13-711-6M	Large	5000



Contact a DCG

Technical Specialist for more information regarding this product.

Call Today!

281-648-1894, Option 1

Sample Vials

Convenient for sample prep or storing samples

<u>Part Number</u>	<u>Size</u>	<u>Per Case</u>
DCG196503	12 mL Clear	200
DCG196503A	12 mL Amber	200
DCG10310-5	<u>Caps</u>	100
DCG15100TS	<u>Septa</u>	100

Septa GC Inlets

Inlet temperature to 400°C Preconditioned

<u>Part Number</u>	<u>Type</u>	<u>Per Pack</u>
DCG27110	11mm Center Guide	50

Graphite Ferrules

For installing capillary columns made of high density, high purity graphite stable to 450°C

<u>Part Number</u>	<u>Type</u>	<u>Per Pack</u>
DCG20250	0.4"	10
DCG21007	0.5"	10

Septum Puller

For removing septa and O-rings also good for removing stuck ferrules or ferrule pieces.

<u>Part Number</u>	DCG20117
--------------------	----------

Ceramic Scoring Wafer

Used to make clean cuts on capillary columns. Important for proper installation of the column in the inlet and detector.

<u>Part Number</u>	DCG20116	5 pk
--------------------	----------	------



DCG Partnership 1, Ltd.

4170A S. Main St. Pearland, TX 77581

281-648-1894, Option 1

www.dcgpartnership.com