


Calibration Standards for ASTM D 4951 or EN 14107: The Elemental Analysis of Phosphorus in Biodiesel

DCG has available a calibration standard for ASTM D 4951 or EN 14107. These test methods are used to determine the Phosphorus content in Biodiesel. These standards are NIST traceable by weight and are analyzed to assure quality. Custom standards are available upon request.

Phosphorus has been shown to inhibit the after treatment process in emission systems. The affect of phosphorus on emission systems is cumulative, so small amounts of phosphorus can be detrimental.

Part #: DCGORGM-P100: Phosphorus (P) at 100 ppm by weight packaged in a 100 mL bottle.



Prices subject to change

4170 A Main Pearland, TX 77581 PH: 281-648-1894 FAX: 281-648-1895
www.dcgpartnership.com
E-mail: customerservice@dcgpartnership.com

Chromatographic Reference Materials