


Calibration Standards for EN 14111: The Analysis of Iodine Value in Biodiesel

DCG has available a calibration standard for **EN 14111**. This test method is used to determine the Iodine Value of Biodiesel fuels. These standards are NIST traceable by weight and are analyzed to assure quality. Custom standards are available upon request.

The Iodine Value is a measure of a biodiesel fuels chemical stability against oxidation.

Part #: DCGEN14111: A set of Standard Iodine Number 50 and Standard Iodine Number 200 packaged in 500 mL bottles. One of each per set.



Prices subject to change

4170 A Main Pearland, TX 77581 PH: 281-648-1894 FAX: 281-648-1895
www.dcgpartnership.com
E-mail: customerservice@dcgpartnership.com

Chromatographic Reference Materials