

DCG Spec-Brief™

SB – 143

Calibration Standards for ASTM D 4951 or EN 14107: The Elemental Analysis of Phosphorus in Biodiesel



DCG has available a calibration Standard for **ASTM D 4951 or EN 14107.** These test methods are used to determine the Phosphorus content in Biodiesel. These standards are N.I.S.T. traceable by weight and are analyzed to assure quality. Custom standards are upon request.

Phosphorus has been shown to inhibit the after treatment process in emission systems. The effect of phosphorus on emission systems is cumulative, so small amounts of phosphorus can be detrimental.

Part #: DCGORGM-P200:

Phosphorus (P) at 100 ppm by weight packaged in a 100 ml. bottle.



DCG PARTNERSHIP 1, LTD.

4170A S. Main Street, Pearland, TX 77581
281-648-1894, Option 1
www.dcgpartnership.com