

DCG Spec-Brief™

SB – 158

Calibration Standards for EN 14214: The Analysis of Iodine Value in Biodiesel



DCG has available a calibration Standard for **EN 14214**. These test methods are used to determine the Iodine Value of Biodiesel fuels. These standards are N.I.S.T. traceable by weight and are analyzed to assure quality. Custom standards are upon request.

The Iodine Value is a measure of a biodiesel fuels chemical stability against oxidation.

Part #: DCGEN14214:

A set of Standard Iodine Number 50 and Standard Iodine Number 200 packaged in 500 ml. bottles. One of each per set.



DCG PARTNERSHIP 1, LTD.

4170A S. Main Street, Pearland, TX 77581
281-648-1894, Option 1
www.dcgpartnership.com